

Validation and Premium Service Plans

Description	Part Number
<p>Comprehensive Twice Yearly IQ/OQ of System with full documentation provided to satisfy regulatory authorities. Comprehensive system test verifies, calibrates and documents all system functions including readings, databases, alarms, dial-outs, outputs, reports, and user programming. Includes One Year Onsite Service Plan with 24/7 Technical Support. Includes approximately 100 page GMP/GLP IQ/OQ validation report, 1000 page software validation test report, calibration of temperature and humidity with traceable equipment, copies of reference equipment calibration certificates and onsite refresher training at PM visit.</p>	<p>VALSERV A VAL/SERV A 3YR Three Year Service Package</p>
<p>Comprehensive Once Yearly IQ/OQ of System with full documentation provided to satisfy regulatory authorities. Comprehensive system test verifies, calibrates and documents all system functions including readings, databases, alarms, dial-outs, outputs, reports, and user programming. Includes One Year Onsite Service Plan with 24/7 Technical Support. Includes approximately 100 page GMP/GLP IQ/OQ validation report, 1000 page software validation test report, calibration of temperature and humidity with traceable equipment, copies of reference equipment calibration certificates and onsite refresher training at PM visit.</p>	<p>VALSERV B VAL/SERV B 3YR Three Year Service Package</p>
<p>Complete IQ/OQ of System. All Inputs Verified and Calibrated, Or Existing Customer Calibrations Rees-Certified.</p>	<p>VALSECOND</p>
<p>Service Contracts: One year service contract, parts and labor, and one year tech service.</p>	<p>SERVICE CONTRACT</p>
<p>One year of technical service phone support.</p>	<p>TECH SERVICE-1</p>

All specifications are subject to change without notice.

One month of technical service phone support.	TECH SERVICE-M
Premium Calibration Services: Five Point Certified Sensors with Certificate (RTD Sensors and High Accuracy Humidity Sensors)	CAL-LOOP
Three Point Calibration (available for those choosing to stay with thermistor probes and RTD's only) Three point calibration service. Field service technicians calibrate sensors using highly accurate dry blocks and probes. Printed certificates of calibrated points will be supplied.	CAL-3POINT
Alarm Limit Testing -More commonly referred to as "Fire and Ice" testing, the Alarm Limit Testing protocol is designed to prove that system will indeed alarm at each sensors set high and low alarm parameters.	FIRE-ICE TEST
Custom onsite calibration points for 3 point calibrations. Required in addition to 3 point calibration part number for each sensor being calibrated, when standard Rees calibration points are not used. Standard Rees calibration points are 937C,25C,-10C,-40C and -70C).	CAL-3POINT-CST
3 Point Onsite Calibration per Temperature Input. Standard Sensor - Range of calibration - 40C to +100C. Includes Calibration Certificate. This option requires that the sensor wire is long enough to be removed from the unit at time of calibration service. If not #CAL-EXT will be required. Minimum of Quantity of 10 required. Note: Requires customer to provide Rees technician remote access to server on a mobile platform such as a laptop. (Note: Requires software Build 1020.15 or higher).	CAL-3POINT-STD
3 Point Onsite Calibration per Temperature Input. Ultra-Low Sensor - Range of calibration - 80C to -40C. Includes Calibration Certificate. This option requires that the sensor wire is long enough to be removed from the unit at time of calibration service If not #CAL-EXT will be required. Minimum of Quantity of 10 required. Note: Requires customer to provide Rees technician remote access to server on a mobile platform such as a laptop. (Note: Requires software Build 1020.15 or higher).	CAL-3POINT-UL
Premium Sensor Calibration Protocols Available	

All specifications are subject to change without notice.

Premium Calibrations

Standard Part Number	Factory Certified Options		Field Certified Options				
	17025- 3 POINT CERT: Calibrated by Rees Scientific certified Presidio Labs and technicians. Includes ISO 17025 certificate.	Add on Cert Available	17025 Single Point: CAL-17025	17025 Multi-Point: CAL-3POINT()	17025 Multi Point: SNAP-CERT-(): Trade out part # with Certified Sensor		
TPT2	TPT2-CERT		CAL-17025	CAL-3POINT-STD	SNAP-CERT-2		
TPT3	TPT3-CERT		CAL-17025	CAL-3POINT-STD	SNAP-CERT-3		
TPT4	TPT4-CERT		CAL-17025	CAL-3POINT-UL	SNAP-CERT-4		
Z3-RTD-STND-1	Z3-RTD-CERT-1		CAL-17025	CAL-3POINT-(STD/UL)			
Z3-RTD-STND-2	Z3-RTD-CERT-2		CAL-17025	CAL-3POINT-(STD/UL)			
Z3-T	Z3-T-CERT		CAL-17025		SNAP-CERT-2		
Z3-T-WF	Z3-T-CERT-WF		CAL-17025		SNAP-CERT-2		
Z3-T-CO2	Z3-T-CO2-CERT				CAL-3POINT-STD CAL-CO2-3PT		
Z3-T-CO2-WF	Z3-T-CO2-WCERT				CAL-3POINT-STD CAL-CO2-3PT		
Z3-T/H	Z3-T/H-CERT *				CAL-CERT-PLUS		SNAP-CERT-TH
Z3-T/H-WF	Z3-T/H-WF-CERT *				CAL-CERT-PLUS		SNAP-CERT-TH
Z3-T/H-REM	Z3-T/H-REM-C *				CAL-CERT-PLUS		SNAP-CERT-TH-R
Z3-T/H-REM-WF	Z3-T/H-RM-WF-C *				CAL-CERT-PLUS		SNAP-CERT-TH-R
HUM-RT1	HUM-RT1-C				Included		SNAP-CERT-RT1
HUM-RT1-REM	HUM-RT1-REM-C	Included			SNAP-CERT-RT1		
SENSMOD-TH	SENSMOD-TH-CRT*	CAL-CERT-PLUS					
SENSMOD-TH-RM	SENSMOD-TH-RMC*	CAL-CERT-PLUS					
SENSMOD-TH-F	SENSMOD-THFCRT*	CAL-CERT-PLUS					

* Includes 3 Pt. 17025 Traceable Certificate for humidity at a single temperature. For 3 point temperature calibration add CAL-CERT-PLUS part number.

Standard Part Number	Factory Certified Options		Field Certified Options		
	17025- 3 POINT CERT: Calibrated by Rees Scientific certified Presidio Labs and technicians. Includes ISO 17025 certificate.	Add on Cert Available	17025 Single Point: CAL-17025	17025 Multi-Point: CAL-3POINT()	17025 Multi Point: SNAP-CERT-(): Trade out part # with Certified Sensor
PSTR-DUCT-()				CAL-DIFF-3PT	
CO2-MON-RM3				CAL-CO2-3PT	
RTD+10+70	RTD+10+70-CRT		CAL-17025	CAL-3POINT-STD	
RTD-0+100	RTD-0+100-CRT		CAL-17025	CAL-3POINT-STD	
RTD-100+30	RTD-100+30-CRT		CAL-17025	CAL-3POINT-UL	
RTD-20+30	RTD-20+30-CRT		CAL-17025	CAL-3POINT-STD	
RTD-200+60	RTD-200+60-CRT		CAL-17025	CAL-3POINT-UL	
RTD-50+50	RTD-50+50-CRT		CAL-17025	CAL-3POINT-STD	
RTD-AMBIENT	RTD-AMBIENTCRT		CAL-17025	CAL-3POINT-STD	

* Includes 3 Pt. 17025 Traceable Certificate for humidity at a single temperature. For 3 point temperature calibration add CAL-CERT-PLUS part number.

International Part Numbers						
Standard Part Number	Factory Certified Options		Field Certified Options			
	17025-3 POINT CERT: Calibrated by Rees Scientific certified Presidio Labs and technicians. Includes ISO 17025 certificate.	Add on Cert Available	17025 Single Point: CAL-17025	17025 Multi-Point: CAL-3POINT()	17025 Multi Point: SNAP-CERT-(): Trade out part # with Certified Sensor	
Z3-RTD-STND1EV	Z3-RTD-CRT1-EV		CAL-17025	CAL-3POINT-(STD/UL)		
Z3-RTD-STND2EV	Z3-RTD-CRT2-EV		CAL-17025	CAL-3POINT-(STD/UL)		
Z3-T-EV	Z3-T-CERT-EV		CAL-17025		SNAP-CERT-2	
Z3-T-WF-EV	Z3-T-WF-EV-C		CAL-17025		SNAP-CERT-2	
Z3-T-CO2-EV	Z3-T-CO2-EV-C				CAL-3POINT-STD CAL-CO2-3PT	
Z3-T-CO2-WF-EV	Z3-T-CO2-WC-EV				CAL-3POINT-STD CAL-CO2-3PT	
Z3-T/H-EV	Z3-T/H-EV-CERT *				CAL-CERT-PLUS	
Z3-T/H-EV-WF	Z3-T/H-EV-WF-C *				CAL-CERT-PLUS	SNAP-CERT-TH
Z3-T/H-REM-EV	Z3-T/H-RM-EV-C *				CAL-CERT-PLUS	SNAP-CERT-TH-R
Z3-T/H-RM-WF-E	Z3-T/H-RM-E-WC *				CAL-CERT-PLUS	SNAP-CERT-TH-R

* Includes 3 Pt. 17025 Traceable Certificate for humidity at a single temperature. For 3 point temperature calibration add CAL-CERT-PLUS part number.

All specifications are subject to change without notice.