



Introducing Rees Scientific's Premier Line of Laboratory Equipment

Lab & Pharmacy Freezers

Our complete line of laboratory and pharmacy freezers provide the reliability and security essential for the protection of sensitive lab samples in the 0°C to -30°C range. Ideal for the demanding requirements of daily lab applications, our range of freezers offer dependability and safety necessary for your daily use.



Lab & Pharmacy Freezers

Our Lab and Pharmacy freezers were designed specifically for use in a laboratory and pharmacy setting. These freezers offer many storage capabilities such as shelves and racks to store samples, chemicals, medications, and other laboratory products. Rees Scientific also offers various options to monitor these freezers for fluctuation in temperature, humidity, differential pressure, etc. New equipment can be easily added to the existing Rees monitoring system. Rees Scientific has set the industry standard for environmental monitoring in pharmacy laboratories for 30 years now. We continue to lead the field, taking cues from our clients and their environmental challenges.

Features:

Laboratory & Pharmacy Freezer

- Digital Temperature Display
- Auto or Manual Defrost Available
- Adjustable Shelves
- CFC-Free Refrigerant and Insulation
- Optional Chart Recorder
- UL Listed

Glass Door Pharmacy Freezer

- Digital Temperature Display
- Auto Defrost
- Adjustable Shelves
- CFC-Free Refrigerant and Insulation
- 3" Casters
- Optional Chart Recorder
- UL Listed

MODEL	CU. FT./LITERS	TEMPERATURE RANGE	INTERIOR DIMENSIONS, INCHES (CM) WxF-BxH	EXTERIOR DIMENSIONS, INCHES (CM) WxF-BxH
LAB & PHARMACY FREEZERS				
DHN40-21SD	21/594	-30°C to -40°C	23.5" x 30" x 59" (59 x 76 x 150)	27.5" x 36" x 81.5" (69 x 90 x 207)
DH29-23SD	23.1/654	-15°C to -30°C	22" x 28" x 60" (56 x 71 x 152)	26" x 33.5" x 84.5" (66 x 85 x 215)
DH29-34SD	34/963	-18°C to -30°C	31" x 28" x 60" (79 x 71 x 152)	35" x 34" x 81.5" (89 x 86 x 207)
DH29-46SD	46.6/1320	-15°C to -30°C	48" x 28" x 60" (122 x 71 x 152)	52" x 33.5" x 84.5" (132 x 85 x 215)
DHN20-52SD	52/1472	-10°C to -25°C	51" x 30" x 59" (129 x 76 x 150)	55" x 35" x 80" (139 x 88 x 202)
DHF29-74SD	74/2095	-15°C to -30°C	74" x 28" x 60" (188 x 71 x 152)	78" x 33.5" x 84.5" (198 x 85 x 215)
DHN20-80SD	80/2265	-10°C to -25°C	78.5" x 30" x 59" (199 x 76 x 150)	82.5" x 35" x 80" (209 x 88 x 202)
DHN20-16MDF	17/484	-10°C to -20°C	27.5" x 21.5" x 52" (69 x 54 x 132)	32" x 30.5" x 64.5" (71 x 77 x 163)
DHN20-24MDSD	24/679	-10°C to -25°C	23.5" x 30" x 59" (60 x 76 x 150)	27.5" x 34.88" x 79.75" (69 x 88 x 202)
DHW20-17MDP	17/480	-15°C to -20°C	27" x 21.5" x 52" (68 x 54 x 132)	32.5" x 29" x 75" (82 x 73 x 190)
DHW20-30MDP	30/849	-15°C to -20°C	32" x 22.5" x 65" (81 x 57 x 165)	35.5" x 35" x 79" (90 x 89 x 200)
DHW20-30	30/849	-15°C to -20°C	32" x 22.5" x 65" (81 x 57 x 165)	35" x 35" x 72" (89 x 89 x 183)
DHW30-30MDP	30/849	-20°C to -30°C	32" x 22.5" x 65" (81 x 57 x 165)	35.5" x 35" x 79" (90 x 89 x 200)
DHW30-30	30/849	-20°C to -30°C	32" x 22.5" x 65" (81 x 57 x 165)	32.5" x 29" x 75" (82 x 73 x 190)

