

Introducing The New Centron SQL

New Build 1014 SQL Software



Rees Scientific's New Software offers Multi-Level Maps providing each department with their own floor maps! See Real time alarms as they occur across the hallway or across the world!

Features Include:

- System wide overview shows ALL nodes locations and status at a glance!
- Click directly to any node, any department or any map on the new overview tree!
- Improved system wide running log allows users to access event overview and comment box simultaneously!

Rees Monitoring System - System Wide Overview

File Alarms Program Window Help

Overview Tree

- Global Overview
 - Node 01-Boston
 - Bld 1 Upper West
 - Bld 1 Upper East
 - Bld 1 Lower West
 - Bld 1 Lower East
 - Node 02-Beijing
 - West Wing Lab
 - Cell Biology
 - Warehouse + Quarantine
 - PD and Virology
 - Northeast Virus Lab
 - Northeast Tissue Lab
 - Southeast Tissue Lab
 - Node 03-Madison
 - Node 04-Madrid
 - Node 05-New Delhi
 - Node 06-Indianapolis
 - Node 07-Austin
 - Node 08-Dallas
 - Node 09-San Diego
 - Node 10-Leeds
 - Node 12-Daytona
 - Node 13-Boise
 - Node 14-Portland
 - Node 15-Honolulu
 - Node 16-Berlin
 - Node 20-Houston
 - Node 21-Sidney Security

System-wide Running Log

Date/Time	In/Out	Node	Condition	Comment
6/ 1/2009 - 14:24	69	Node 06-Indianapolis	Program Change	Yes
6/ 1/2009 - 14:25	69	Node 06-Indianapolis	Program Change	Yes
6/ 1/2009 - 14:25	69	Node 06-Indianapolis	-67.6 In Alarm	No
6/ 1/2009 - 14:25	69	Node 06-Indianapolis	-67.6 Phoning	No
6/ 1/2009 - 14:33	69	Node 06-Indianapolis	Disable Alr 1.00hr	No
6/ 1/2009 - 14:35	69	Node 06-Indianapolis	Alarm Reset	No
6/ 1/2009 - 14:35	69	Node 06-Indianapolis	-67.6 In Alarm	No
6/ 1/2009 - 14:35	69	Node 06-Indianapolis	-67.6 Phoning	No
6/ 1/2009 - 14:44	80	Node 02-Beijing	Program Change	Yes
6/ 1/2009 - 15:05	107	Node 02-Beijing	Program Change	Yes
6/ 1/2009 - 15:05	107	Node 02-Beijing	Program Change	Yes

System Wide Overview

Node 2: Node 02-Beijing

Node 4: Node 04-Madrid

Ready

ONLINE CAP NUM

Rees Scientific

1007 Whitehead Road Ext. Trenton, NJ 08638 P(609)530-1055 F(609)671-2751

Email: Sales@reesscientific.com www.reesscientific.com

Introducing The New Centron SQL



Departmental Reviews and Approval - This new optional module will allow each department to perform review and approval of their own inputs. Previously, Supervisory permissions only allowed review and approval of ALL the inputs. This new feature will allow Organizations to delegate this responsibility to each department. Additionally, the system will no longer require a Supervisor Security clearance to perform the review and approval process, allowing Managers to delegate this responsibility to others within the department.

Security Code and Telephone Programming - Previously, the software required that ALL Nodes be open when making changes to Security Codes or Telephone programming. This often caused an issue when a user did not have access to Nodes in other departments or buildings and therefore, could not perform the changes. It also required additional time for the system to open all the Nodes. The new software will allow users to make these changes to only their Node or Nodes without needing access to other Nodes on the system.

