



Introducing Rees Scientific's Premier Line of Laboratory Equipment

Undercounter Freezers / Refrigerators

Our complete line of undercounter refrigerators and freezers provide the reliability and security essential for the protection of your sensitive commodities. Manufactured with durable commercial grade construction, our Undercounter Freezers and Refrigerators cool quickly and efficiently to avoid temperature fluctuations which can result in failure of expensive research, misleading test results, or weakened medicines, vaccines, and reagents.



Undercounter Refrigerators / Freezers

These general purpose undercounter models can be built into casework and under lab benches, offering accessibility and superior storage opportunities. A wide array of options and accessories are available, including lockable storage inserts, stainless steel cabinet exterior, and glass doors. Rees Scientific also offers various options to monitor these refrigerators and freezers for fluctuation in temperature, humidity, differential pressure, etc. New equipment can be easily added to the existing Rees monitoring system. Rees Scientific has set the industry standard for environmental monitoring for 30 years now. We continue to lead the field, taking cues from our clients and their environmental challenges.

Features:

Undercounter Refrigerators / Freezers

- Durable Construction
- Manual Defrost and Auto Defrost
- Epoxy Coated Shelving
- Air Cooled Condenser
- Temperature Fluctuation Reduction
- CFC-Free Refrigeration System
- Undercounter Storage Capabilities
- Lockable Storage Inserts
- Stainless Steel Versions
- UL and CUL Listed

MODEL	CU. FT./LITERS	TEMPERATURE RANGE	INTERIOR DIMENSIONS, INCHES (CM) WxF-BxH	EXTERIOR DIMENSIONS, INCHES (CM) WxF-BxH
UNDERCOUNTER REFRIGERATORS & FREEZERS				
MV25-4UCF	4.5/127	-17°C to -25°C	(FREEZER WITH ALARM)	24" x 24" x 33.5"
MV30-4UCF	4.5/127	-20°C to -30°C	(FREEZER WITH ALARM)	24" x 24" x 33.5"
MV23-4UCF	4.5/127	-18°C to -23°C	(FREEZER)	24" x 24" x 34"
MV12-6UCF	6.1/173	-7°C to -12°C	(FREEZER)	24" x 24" x 34"
MV20-1.5UCF	1.5/42	-20°C	(FREEZER)	15" x 22.5" x 24.5"
MV4-6UCRDA	6.1/173	1°C to 7°C	(REFRIGERATOR WITH ALARM)	24" x 24" x 34"
MV4-6UCRGDA	6.1/173	1°C to 7°C	(REFRIGERATOR WITH GLASS DOOR)	24" x 24" x 34"
MV4-6UCRGDDA	6.1/173	1°C to 7°C	(REFRIGERATOR WITH GLASS DOOR)	24" x 24" x 34"
MV4-6UCRGDA-ADA	6.1/173	1°C to 7°C	(REFRIGERATOR WITH GLASS DOOR AND ALARM)	24" x 24" x 34"
MV4-8UCRDA	8/226	1°C to 7°C	(REFRIGERATOR WITH GLASS DOOR)	24" x 24" x 34"
MV18-6UCRF	5.3/150	1°C to 7°C	(REFRIGERATOR / FREEZER)	24" x 24" x 34"
	0.7/23	-9°C to -18°C		
MV23-4UCXF	4.5/127	-15°C to -18°C	(EXPLOSION PROOF FREEZER)	24" x 24" x 34"
MV23-4UCFMSF	4.5/127	-15°C to -18°C	(FLAMMABLE STORAGE FREEZER)	24" x 24" x 34"
MV4-6UCXR	6.1/173	1°C to 7°C	(EXPLOSION PROOF REFRIGERATOR)	24" x 24" x 34"
MV4-6UCFMSR	6.1/173	1°C to 7°C	(FLAMMABLE STORAGE REFRIGERATOR)	24" x 24" x 34"
DHN4-5UCR	5.1/144	2°C to 8°C	(REFRIGERATOR)	27" x 25" x 34"
DHN4-5MDUCR	5.1/144	2°C to 8°C	(REFRIGERATOR)	27" x 25" x 34"
DHN20-4UCF	5.1/144	-10°C to -20°C	(FREEZER - AUTO DEFROST)	27" x 25" x 34"
DHN20-4MDUCF	5.1/144	-10°C to -20°C	(FREEZER)	27" x 25" x 34"
DHN30-5UCF	5.1/144	-30°C	(FREEZER - AUTO DEFROST)	27" x 25" x 34"

