

Rees Scientific is proud to introduce our New 5 GHz WiFi system



Our **5 GHz Wi-Fi** offers the same great reliability, security and flexibility as the Z3 Wireless System. It's designed to reduce the overall installation process and is great for remote sites. It offers enhanced customer convenience and works with the full line of Rees Scientific sensors. **Our 5 GHz radios are FIPS (Federal Information Processing Standard) certified.**



The 5 GHz spectrum offers improved radio communication and less interference



Has the ability to operate on 2.4 and 5 GHz

Min 2.0°C
Max 8.0°C

Modules with displays show the daily Min and Max in compliance with vaccine storage requirements



Provides a significantly longer and more reliable battery life under most Wi-Fi conditions



It offers a higher level of security with a FIPS 140-2 Level 1 validated radio



Offers dependable and consistent data buffering



Has a superior ability to stay online even in the most difficult of conditions



Enterprise Wi-Fi capable option (PEAP-MS-CHAPv2)





More security and flexibility with the same great service!



MONITORING



PROTECT VALUABLE ASSETS



COMPLIANCE



NOTIFICATION



SECURED WIFI

