

# THE ENTERPRISE CENTRON SOLUTION



## About

The Enterprise Centron offers a more powerful, reliable and scalable centralized monitoring system. It is expandable up to 512 inputs and employs RAID technology that duplicates data across multiple drivers to provide full redundancy. Two disks each store exactly the same data, at the same time, all the time. Your vital data is protected by dual disks. Strengthen your facilities centralized monitoring system with the Enterprise Centron.

## Features

- \* More powerful & reliable
- \* Multiple dialers
- \* Expandable up to 512 inputs
- \* Simultaneous networking for up to 30 users
- \* RAID dual level data redundancy
- \* Access via web browser
- \* Available in standalone tower or rack mount



MONITORING



PROTECT VALUABLE ASSETS



COMPLIANCE



REGIONALIZED SERVICE



NOTIFICATION



GRAPH OF EVENTS

