



Automated Temperature Monitoring



Joint Commission

You know Joint Commission accreditation is important for improvement in the quality of the staff and quantity of patients and contracts. Seeking accreditation can be a daunting task.

Joint Commission Standards PC.17.10.6-8 require organizations that store or issue tissue, such as pharmacies and blood and tissue banks, to *maintain continuous temperature monitoring for storage refrigerators and freezers as well as maintain daily records to show that tissues were stored at the required temperatures.* It is also a requirement that *storage equipment has functional alarms and emergency back up.*

Rees Scientific's CentronEMS can make the process easier through:

- Alarm Notification**
- Automated Monitoring**
- Compliance**
- Protecting your valuable data**
- Reporting**
- Record Keeping**

In addition to Joint Commission, Rees Scientific's CentronEMS meets today's toughest regulatory requirements including:

- AABB**
- CAP**
- FDA 21 CFR 11**
- GxP**
- HACCP**
- USP 797**



The next generation of Rees Scientific Wireless Environmental Monitoring brings you increased reliability and clarity!

V.2 offers:

- Data buffering at Transmitter Level**
- Battery Alarm**
- Reliable Two-way Communication**
- Optional Local sonic Alarm**

Protecting your critical equipment throughout the hospital

- Blood Banks**
- IVF Labs**
- Microbiology**
- Nutrition**
- Pathology**
- Pharmacy**
- Tissue Banks**



1007 Whitehead Road Ext.
Trenton, NJ 08638 USA
609 530.1055 800 371.3141
www.reesscientific.com
sales@reesscientific.com

PRINT2543.0901

wireless
Rees Scientific V.2
Environmental Monitoring Systems
Integrated Access Control
Automated Watering Solutions

Multiple departments. One answer.

Pharmacy, Blood and Tissue Banks, Nutrition, Core Labs

You all have critical equipment (refrigerators, freezers, incubators and more) that need to be monitored and maintained. Each department's staff spends a considerable portion of their day recording temperature, humidity, pressure and flow readings. Automating environmental monitoring increases accuracy, provides early warnings to address problems before they occur, maintains JCAHO and other regulatory compliance, and allows staff to focus on other critical tasks.

Monitor your equipment regardless of its location. A single Centron system can monitor numerous departments within your facility whether they are on a different floor, in a different building, or even in a different town.

Share costs but not control.

Each input can be assigned to a specific department and viewed only by the authorized users. Allow each department to handle its own alarms, generate its own reports and use its own SOPs.

Module-based design provides flexibility and expandability.

Wireless sensors make retrofitting existing facilities easy and economical and can be used in conjunction with wired sensors.

Address problems before they occur.

With the latest release of our Centron EMS software, you will find new features including a pre-alarm alert system with sonic, relay and email notification providing early warnings that allow you to resolve a potential problem before anything is compromised.

REES SCIENTIFIC.
Can you afford not to?

