



## Introducing Rees Scientific's Premier Line of Laboratory Equipment

### *Lab & Pharmacy Refrigerators*

Our complete line of laboratory and pharmacy refrigerators provides the reliability and security essential for the protection of sensitive lab samples in the +4°C range. Ideal for the demanding requirements of daily lab applications, our range of refrigerators offer dependability and safety necessary for your daily use.



# Lab & Pharmacy Refrigerators

Our Lab and Pharmacy Refrigerators were designed specifically for use in a laboratory and pharmacy setting. These refrigerators offer many storage capabilities such as shelves and racks to store samples, chemicals, medications, and other laboratory products. Rees Scientific also offers various options to monitor these freezers for fluctuation in temperature, humidity, differential pressure, etc. New equipment can be easily added to the existing Rees monitoring system. Rees Scientific has set the industry standard for environmental monitoring in pharmacy laboratories for 30 years now. We continue to lead the field, taking cues from our clients and their environmental challenges.

## Features:

### Laboratory & Pharmacy Refrigerators

- LCD Temperature Readout
- High & Low Alarm with Dry Contacts
- Casters
- Double-Paned Tempered Glass Doors
- Adjustable Shelves
- Adjustable Thermostat
- Cycle Defrost
- CFC-Free Refrigeration System

MODEL	CU. FT./LITERS	TEMPERATURE RANGE	INTERIOR DIMENSIONS, INCHES (CM) WxF-BxH	EXTERIOR DIMENSIONS, INCHES (CM) WxF-BxH
<b>LAB &amp; PHARMACY REFRIGERATORS</b>				
DH4-10GD	10/283	1°C to 12°C	21" x 20.75" x 37.75" (52 x 53 x 96)	24" x 24" x 65" (60 x 60 x 165)
DH4-12GD	12/339	1°C to 12°C	21" x 20.75" x 45.75" (52 x 53 x 116)	24" x 24" x 73" (60 x 60 x 185)
DH4-15SGD	15/424	1°C to 12°C	33" x 21" x 40.5" (83 x 53 x 102)	36" x 24" x 65" (91 x 60 x 165)
DH4-17SGD	17/481	1°C to 12°C	33" x 21" x 46.5" (83 x 53 x 118)	36" x 24" x 71" (91 x 60 x 180)
DH4-23GD	23/651	1°C to 12°C	24" x 28.5" x 54" (60 x 72 x 137)	27" x 32" x 81.5" (68 x 81 x 207)
DH4-27GD	23/764	1°C to 12°C	27" x 28.5" x 55" (69 x 72 x 140)	30" x 32" x 81.5" (76 x 81 x 198)
DH4-27SD	27/764	1°C to 12°C	27" x 28.5" x 55" (69 x 72 x 140)	30" x 32" x 81.5" (76 x 81 x 198)
DH4-38SGD	38/1076	1°C to 12°C	40.5" x 25.5" x 55" (102 x 65 x 140)	43.5" x 32" x 81.5" (110 x 81 x 198)
DH4-45SGD	45/1286	1°C to 12°C	49" x 25.5" x 55" (125 x 65 x 140)	52" x 32" x 81.5" (132 x 81 x 198)
DH4-49GD	49/1387	1°C to 12°C	49" x 28.5" x 55" (125 x 72 x 140)	52" x 32" x 81.5" (132 x 81 x 198)
DH4-49SD	49/1387	1°C to 12°C	49" x 28.5" x 55" (128 x 72 x 140)	52" x 32" x 81.5" (132 x 81 x 198)
DH4-66SGD	66/1868	1°C to 12°C	72" x 27" x 55" (183 x 68 x 140)	75" x 32" x 81.5" (132 x 81 x 207)
DH4-74GD	74/2095	1°C to 12°C	72" x 28.5" x 55" (183 x 72 x 140)	75" x 32" x 81.5" (191 x 81 x 198)
DH4-74SD	74/2095	1°C to 12°C	72" x 28.5" x 55" (183 x 72 x 140)	75" x 32" x 81.5" (191 x 81 x 198)
DHN4-80GD-S	80/2265	2°C to 10°C	78.5" x 30" x 59" (199 x 76 x 150)	82.5" x 35.5" x 81.5" (209 x 90 x 207)
DHN4-80SD-S	80/2265	2°C to 10°C	78.5" x 30" x 59" (199 x 76 x 150)	82.5" x 35.5" x 81.5" (209 x 90 x 207)

