



MONITORING



PROTECT CRITICAL EQUIPMENT



COMPLIANCE



NOTIFICATION




GRAPH OF EVENTS

## Continuous Temperature Monitoring

### New RTD Temperature Sensor

The same company that has provided protection and peace of mind for over 30 years is introducing their 3 wire precise and accurate Resistance Temperature Detector sensor. Our new RTD sensor comes equipped with a high accuracy wide range probe and one sensor covers the full laboratory temperature range.

- \* Equipped with a PT100 platinum RTD probe
- \* Operates over a -200° C to +100° C range
- \* Offers a .07 degree resolution
- \* Single or 2 input models
- \* Standard version has typical 0.4 degree accuracy over span. Certified version has 0.1 degree accuracy over span (powered).
- \* Certified version comes with a certificate showing 0.1 degree accuracy at -196, -70, -40, -20, 0, 20, and 50 degrees C.
- \* RTDs offer greater stability and accuracy than thermocouples or thermistors.





MONITORING



PROTECT CRITICAL EQUIPMENT



COMPLIANCE



NOTIFICATION



GRAPH OF EVENTS

## Continuous Temperature Monitoring

### How it works

