



MONITORING



PROTECT VALUABLE ASSETS



COMPLIANCE



NOTIFICATION




GRAPH OF EVENTS

## Continuous Temperature Monitoring

### RTD Temperature Sensor

The same company that has provided protection and peace of mind for over 35 years is introducing their 3 wire precise and accurate Resistance Temperature Detector sensor. Our new RTD sensor comes equipped with a high accuracy wide range probe and one sensor covers the full laboratory temperature range.

- \* Equipped with a PT100 platinum RTD probe.
- \* Operates over a -200° C to +100° C range.
- \* Offers a .07 degree resolution.
- \* Single or 2 input models.
- \* Available in Z3 Wireless and WiFi.
- \* Standard version has typical 0.4 degree accuracy over span. Certified version is calibrated to 0.1 degree accuracy at -196, -50, 0, and 50 degrees C. A certificate is provided.
- \* RTDs offer greater stability and accuracy than thermocouples or thermistors.





MONITORING



PROTECT VALUABLE ASSETS



COMPLIANCE



NOTIFICATION



GRAPH OF EVENTS

## Continuous Temperature Monitoring

### How it works

