




# REES DIALER OPTIONS



Back View

# CELLULAR



Rees Scientific has improved and expanded their capabilities and offer more options for their alarm dialing.

## VOIP Dialer

- \* For use with virtualized Centron Presidio installation
- \* Constant connection to box
- \* Requires an analog phone line for operation
- \* Highest level of security for mission critical assets and events
- \* Optional 1 or 4 line dialer

## VOIP Software

- \* Viable option for sites that don't own or want POTS line and TAPI modem
- \* Ability to leverage modern network infrastructure
- \* Less hardware required (eliminates need for modem, telephone line and internal connection of modem)

## VOIP Backup Dialer

- \* Monitors phone line for connectivity
- \* Additionally, if the main phone line service is interrupted, the VOIP backup dialer will take over alarm duty
- \* Another safeguard ensuring your critical equipment is safe and protected

## Cellular Dialer

- \* External Cellular Dialer to connect to customer provided PC
- \* Performs telephone alarm dial-out for the Rees Scientific Virtual Node and requires 1 external USB port
- \* Synthesized voice dial out capable of working with RSC customizable voice tags
- \* Comes with prepaid GSM SIM Card, but ongoing cellular service through a Cellular GSM telephone provider is required
- \* Can be main dialer or back up