

The new E-CENTRON!
A WIDE ARRAY of FEATURES
in ONE CONVENIENT PACKAGE!

The new E-Centron will combine several exceptional Rees Scientific features into one convenient and economical package.

E-Centron will combine Centron hardware, Web Access, internal CD-RW drive, Autobackup software, 4-hour battery backup, printer, monitor, keyboard and mouse all delivered in a Tower Chassis.



E-Centron will be licensed to accept 16 inputs, wired and/or wireless, with the option of expanding in increments of 16 inputs or licenses. The E-Centron will allow an unlimited number of users to access the system through an off-the-shelf browser via your corporate intranet or the World Wide Web. Users can also interface with the system at the main E-Centron Tower.

Rees Scientific