

Satellite Tower Remote Site Monitor

The **Satellite Tower** is invaluable when you have to oversee multiple locations with only a few input points at each site.

The Satellite Tower will allow you to monitor up to eight inputs and will work with your existing Centron Windows® based software; the Satellite will simply be an additional node on your system. The Satellite Tower has all the capabilities of the Centron system with the exception of output controls, and it will operate over your existing TCP/IP WAN, allowing you to access information from multiple sites.



You have real-time readings of each unit at each site, as well as alarm notification. You also have the ability to print out reports directly at the main facility.



- **The Satellite is a fully validated 12 bit system.**
- **It works with all Rees Scientific pluggable sensors.**
- **It has off-site dialing capabilities.**
- **It operates over your existing TCP/IP network.**
- **Completely stand-alone w/4 hour battery backup.**